



鼻梁型 (I)									单位 (CM)	
		型号	A	B	C	D	E	F	G	用途
 <p>“D”表示从外部边缘至内侧凹陷部位（虚线标示部分）的尺寸</p>		K9010VS-1	4.6	2.5	0.8	0.4	0.6	0.5	0.15	隆鼻
		K9010F-1	4.6	2.5	0.8	0.4	0.6	0.5	0.15	隆鼻
		K9010VS-2	4.8	2.5	0.8	0.4	0.6	0.5	0.15	隆鼻
		K9010F-2	4.8	2.5	0.8	0.4	0.6	0.5	0.15	隆鼻
		K9010VS-3	5	2.5	0.9	0.5	0.6	0.5	0.15	隆鼻
		K9010F-3	5	2.5	0.9	0.5	0.6	0.5	0.15	隆鼻
		K9010VS-4	5.2	2.5	0.9	0.5	0.6	0.5	0.15	隆鼻
		K9010F-4	5.2	2.5	0.9	0.5	0.6	0.5	0.15	隆鼻
		K9010VS-5	5.5	3	1	0.6	0.6	0.5	0.15	隆鼻
		K9010F-5	5.5	3	1	0.6	0.6	0.5	0.15	隆鼻

OF 表示所使用的原料牌号 MED-4840, 硬度: 硬 不透明

OVS 表示所使用的原料牌号 MED-4234, 硬度: 非常软 不透明

鼻梁型 (II)									单位 (CM)
		型号	A	B	C	D	E		用途
		K9020VS-1	406	1.3	0.8	0.4	0.6		隆鼻
		K9020F-1	406	1.3	0.8	0.4	0.6		隆鼻
		K9020VS-2	407	1.3	0.8	0.4	0.6		隆鼻
		K9020F-2	407	1.3	0.8	0.4	0.6		隆鼻
		K9020VS-3	409	1.3	0.9	0.4	0.6		隆鼻
		K9020F-3	409	1.3	0.9	0.4	0.6		隆鼻
		K9020VS-4	5.1	1.3	0.9	0.5	0.6		隆鼻
		K9020F-4	5.1	1.3	0.9	0.5	0.6		隆鼻
		K9020VS-5	5.3	1.3	0.9	0.5	0.6		隆鼻
		K9020F-5	5.3	1.3	0.9	0.5	0.6		隆鼻
		K9020VS-6	5.5	1.3	1	0.6	0.6		隆鼻
		K9020F-6	5.5	1.3	1	0.6	0.6		隆鼻

鼻梁型(III)									单位 (CM)	
		型号	A	B	C	D	E	F	G	用途
		K9130VS-1	4.6	2.8	0.85	0.55	0.33	0.45	0.5	隆鼻
		K9130F-1	4.6	2.8	0.85	0.55	0.33	0.45	0.5	隆鼻
		K9130VS-2	4.7	2.8	0.92	0.6	0.45	0.4	0.65	隆鼻
		K9130F-2	4.7	2.8	0.92	0.6	0.45	0.4	0.65	隆鼻
		K9130VS-3	5	3	1	0.65	0.5	0.4	0.7	隆鼻
		K9130F-3	5	3	1	0.65	0.5	0.4	0.7	隆鼻

鼻脊型(I)							单位 (CM)
		型号	A	B	C	D	用途
		K9030VS-1	4.6	0.8	0.4	0.6	隆鼻
		K9030F-1	4.6	0.8	0.4	0.6	隆鼻
		K9030VS-2	4.7	0.8	0.4	0.6	隆鼻
		K9030F-2	4.7	0.8	0.4	0.6	隆鼻
		K9030VS-3	4.9	0.9	0.4	0.6	隆鼻
		K9030F-3	4.9	0.9	0.4	0.6	隆鼻
		K9030VS-4	5.1	0.9	0.5	0.6	隆鼻
		K9030F-4	5.1	0.9	0.5	0.6	隆鼻
		K9030VS-5	5.3	0.9	0.5	0.6	隆鼻
		K9030F-5	5.3	0.9	0.5	0.6	隆鼻
		K9030VS-6	5.5	1	0.6	0.7	隆鼻
		K9030F-6	5.5	1	0.6	0.7	隆鼻