



鼻窦假体 (I)



型号	A	B	C	用途
K6010VS-1	1.7	1.4	0.3	鼻唇沟手术
K6010VS-2	2	1.7	0.3	鼻唇沟手术
K6010VS-3	2.3	1.9	0.4	鼻唇沟手术
K6010VS-4	2.6	2	0.5	鼻唇沟手术
K6010VS-5	2.7	2.2	0.5	鼻唇沟手术
K6010VS-6	2.8	2.3	0.6	鼻唇沟手术

鼻窦假体 (II)



型号	A	B	C	用途
K6020VS-1	1.7	1.4	0.3	鼻唇沟手术
K6020VS-2	1.9	1.7	0.3	鼻唇沟手术
K6020VS-3	2.2	2	0.4	鼻唇沟手术
K6020VS-4	2.5	2.1	0.5	鼻唇沟手术
K6020VS-5	2.8	2.5	0.5	鼻唇沟手术
K6020VS-6	3.2	2.8	0.6	鼻唇沟手术

鼻窦假体 (III)



型号	A	B	C	用途
K6030VS-1	3	1.6	0.4	鼻唇沟手术
K6030VS-2	3.4	1.8	0.5	鼻唇沟手术
K6030VS-3	3.6	1.9	0.5	鼻唇沟手术
K6030VS-4	3.8	2	0.6	鼻唇沟手术
K6030VS-5	4	2.1	0.6	鼻唇沟手术
K6030VS-6	4.2	2.2	0.7	鼻唇沟手术